

## Delivering independent, science and outcome-based pain relief and muscle rehabilitation

RS Medical seeks to improve the quality of life for patients by providing affordable, independent, specialty pain relief and rehabilitation treatment options with timely, effective service. Since 1990, the company has helped more than 1.5 million patients suffering from painful muscle conditions. Together, we can improve patient outcomes, advance the patient experience and reduce the cost of care.



# Electrotherapy

## RS-4i<sup>®</sup> Plus Sequential Stimulator with Intersperse<sup>®</sup> HCPCS E1399

### Multiple Benefits in Each Treatment

The RS-4i Plus combines high frequency interferential current for pain relief and muscle stimulation for rehabilitation in a patent-pending treatment program known as Intersperse<sup>®</sup>, which delivers both therapies in a single, customizable treatment session that is effective, comfortable and convenient. The RS-4i Plus is a non-addictive alternative to common drug therapies, including opioids.



**Before Treatment:**  
Painful muscle condition



**Segment 1:**  
5 minutes of pain relief



**Segment 2:**  
5 minutes of muscle therapy



**Segment 3:**  
25 minutes of Intersperse<sup>®</sup>, which combines pain relief and muscle therapy

### More Advanced Therapy Than TENS

The sequential stimulation offered in an RS-4i Plus goes far beyond standard TENS therapy. It has four channels that provide coverage for large treatment areas like the upper and the lower back. At 5000 Hz in interferential output intensity, pain relief from the RS-4i Plus may last hours after the treatment, meaning patients use the device just once or twice daily. The inclusion of NMES stimulation provides science-based benefits to improve muscle health that surpass the pain masking of TENS.

### Patients Report Fast-Acting Relief

The RS-4i Plus has a documented record of success with VA patients. Because it is designed for use in a home setting and is easily transported, Veterans can fit treatments into their daily schedule when convenient. Empowering patients to be a part of their care can lead to quicker recovery and more positive outcomes. After one month of use, 84 percent of RS-4i Plus patients report a medically meaningful reduction in pain and one-third report activity improvement of greater than 50 percent.<sup>1</sup>

<sup>1</sup>2018 RS Medical, pre-publication data. Patient reported response to RS-4i Plus treatment; analysis of greater than 5,600 patient assessments from U.S. Veteran treated with the RS-4i Plus and followed for one month.

	TENS	RS-4i <sup>®</sup> Plus
# Channels	2	4
<b>Average Daily Treatment</b>		
Dose	Continuous use	35-min treatment
Typical use per day	1-12 hours	1-2 times/day
<b>Indications for Use</b>		
Relieve acute pain	X	X
Relieve chronic pain	X	X
Retard disuse atrophy		X
Reduce muscle spasm		X
Increase local blood circulation		X
Increase range of motion		X
Re-educate muscle		X


SIN	Catalog Number	Product Description	Award Price
A-72	950-0036-00	RS-4i <sup>®</sup> Plus Sequential Stimulator	\$1402.94 ea.

A total of eight (8) electrode pads are included for initial use. Humidity, skin preparation and storage will age electrodes. Monthly resupply is 2 packs of electrode pads.





# Wearable Garments



## RS-LBG™ Low Back Garment **HCPCS E0731**

Accessory garment consistently and accurately applies all eight electrodes on difficult to reach low back locations as easily as putting on a waist belt. Design adjusts in size and utilizes light, comfortable fabrics to fit under clothing for discreet use. The Low Back Garment is made for the RS-4i Plus.

SIN	Catalog Number	Product Description	FSS Award Price
A-72	959-0036-00	RS-LBG® Low Back Garment	\$192.54 ea.

A total of eight (8) electrode pads are included for initial use. Humidity, skin preparation and storage will age electrodes. Monthly resupply is 2 packs (4 pads per pack) of electrode pads.



An example of a full-back 4-channel bipolar electrode placement

## RS-FBG™ Full Back Garment **HCPCS E0731**

Accessory garment consistently and accurately applies all eight electrodes on difficult to reach upper and middle back locations as easily as putting on a vest. Design adjusts in size and utilizes light, comfortable fabrics to fit under clothing for discreet use. The Full Back Garment is made for the RS-4i Plus.

SIN	Catalog Number	Product Description	FSS Award Price
A-72	959-0032-00	RS-FBG® Full Back Garment	\$192.54 ea.

A total of eight (8) electrode pads are included for initial use. Humidity, skin preparation and storage will age electrodes. Monthly resupply is 2 packs (4 pads per pack) of electrode pads.



## Electrode Pads **HCPCS A4556**

Self-adhesive, reusable electrodes of various size and shapes. Designed for use with RS Medical electrotherapy devices and conductive garments to ensure consistent, safe, comfortable electrotherapy treatments. Each pack contains 4 electrodes.

### Standard Electrodes

SIN	Catalog Number	Product Description	FSS Award Price
A-72	180-0011-00	2" Round Electrodes (pack of 4)	\$11.78
A-72	180-0009-00	2.75" Round Electrodes (pack of 4)	\$11.78
A-72	180-0010-00	2" Square Electrodes (pack of 4)	\$11.78
A-72	180-0007-00	2" Skin-Friendly Electrodes (pack of 4)	\$11.78

### Wearable Garment Electrodes

SIN	Catalog Number	Product Description	FSS Award Price
A-72	180-0040-00	2" Square VELCRO® Top Electrodes (pack of 4) <b>For Use With Garments Purchased AFTER 10-1-2018</b>	\$11.78
A-72	180-0005-00	2" Square Double-Sided Garment Electrodes (pack of 4) <b>For Use With Garments Purchased BEFORE 9-30-2018</b>	\$11.78

# Cervical Traction



## Pronex® Inline Cervical Traction **HCPCS E0855**

Unique cervical traction device that supports the natural curvature of the cervical spine, gently stretches soft tissue and provides an even distraction in the anterior and posterior cervical discs for the management of chronic cervical pain.



SIN	Catalog Number	Product Description	FSS Award Price
A-70	802-0011-00	Pronex / Pneumatic Cervical Traction (REG: 14" - 16")	\$226.12 ea.
A-70	802-0012-00	Pronex / Pneumatic Cervical Traction (L: 16" - 18")	\$225.76 ea.
A-70	802-0013-00	Pronex / Pneumatic Cervical Traction (W: 18" + )	\$225.12 ea.



## ComforTrac™ Cervical Traction **HCPCS E0849**

Comfortable, easy-to-use device delivers precise traction with a gauged hand pump and easily adjustable slide stand for targeted relief of chronic cervical pain.

Open Market Price: **\$340.00 ea.**

Add \$25.00 ea. for duffel bag

## VA/FSS Award Details

International Rehabilitative Sciences, Inc. dba RS Medical

Contract # 36F79720D0003

Award Dates: 11/01/2019 - 10/31/2024

65 II A Medical Equipment & Supplies

FSC Class: 6515

NAICS Code: 334510

DUNS # 625133228

Tax ID # 91-1490213

Structure:

Commercial

FOB:

Destination

Payment Address:

PO Box 872650

Vancouver, WA 98684 USA

Purchase Cards:



Fax: 888-496-0121

Email: [va@rsmedical.com](mailto:va@rsmedical.com)

Questions: Speak with an RS Medical VA Specialist at 866-832-4408

**RS Medical** 

14001 SE First St | Vancouver, WA 98684-3509

[www.rsmedical.com](http://www.rsmedical.com)

<sup>1</sup>2012 RS Medical, pre-publication data. Patient reported response to electrotherapy treatment; Analysis of 445,000 patient surveys. Subset analysis of 1,300 person compliant patient cohort with lumbago, treated with the RS-4i and followed for over 6 months. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician. See individual product Instructions For Use for full prescribing information. All products and brand names are trademarks or registered trademarks of their respective companies or organizations.

©2022 RS Medical. All rights reserved. 331-0091-00 Rev. 1 09/22

