

The ComforTrac™ and Pronex® cervical traction devices are patient-controlled, at-home therapies intended to treat musculoskeletal or neurologic impairment of the cervical spine to relieve pain, relax muscle spasms and decompress spinal structures.

ComforTrac™

- ComforTrac includes a patented dual-action hand pump with gauge and quick-release button to make it easy to adjust traction force
- Provides up to 50 lb. of force
- Easily adjustable slide stand with 10, 15 and 20 degree angles
- Contoured memory foam headrest
- Independently adjustable neck wedges provide a custom fit

Pronex®

- Pronex provides even distraction in the anterior and posterior cervical discs
- Supports and reinforces the natural curvature of the spine¹
- Reduces headaches from excessive pressure on the occiput
- Air-inflated bellows provides up to 35 lb. of force
- Flexion wedge allows for 15 and 22 degree angles
- No strain on the jaw's temporomandibular joint (TMJ)

Pronex		
Size	Measurement	
Regular	14-16"	Fits most women of typical height, weight and neck size
Large	16-18"	Fits most men of typical height, weight, and neck size
Wide	18" +	Fits individuals with large shoulders and neck or a short neck



RS Medical  More practice.
Less pain.

14001 SE First Street | Vancouver, WA 98684 | Toll Free: 1-855-490-PAIN (7246) | www.rsmedical.com